Manufacturer  
Fabricante  
Hersteller  
Fabricant  
Fabbricante  
Valmistaja  
製造商  
製造元  
Изготовитель  
Producător  
Tillverkare  
Производитель  
Proizvođač  
Výrobce  
Producent  
Fabrikant  
Κατασκευαστής  
Gyártó  
Produsent  
Producent  
Fabricante  
Proizvodáč  
Výrobca  
Izdelovalec  
Üretici  
الجهة المصنعة

Date of Manufacture  
Fecha de fabricación  
Herstellungsdatum  
Date de fabrication  
Data di produzione  
Valmistuspäivämäärä  
製造日期  
製造年月日  
Дата изготовления  
Data fabricației  
Tillverkningsdatum  
Дата на производство  
Datum proizvodnje  
Datum výroby  
Fremstillingsdato  
Productiedatum  
Ημερομηνία κατασκευής  
A gyártás dátuma  
Produksjonsdato  
Data produkcji  
Data de fabrico  
Datum proizvodnje  
Dátum výroby  
Datum izdelave  
Üretim Tarihi  
تاريخ التصنيع

CE Certification Mark  
Marca de certificación CE  
CE Zertifizierungsmarke  
Marquage CE  
Marcatura CE  
CE-sertifiointimerkki  
CE 认证标记  
CE 認証マーク  
Знак сертификации CE  
Marcă de certificare CE  
CE-märkning  
CE сертификационная маркировка  
Certifikacijska oznaka CE  
Označení CE certifikace  
CE-mærke  
CE-markering  
Σήμα πιστοποίησης CE  
CE tanúsítvány-jelölés  
CE-sertifiseringsmerke  
Znak CE  
Marca de certificação CE  
CE oznaka  
Certifikačná známka CE  
Oznaka CE  
CE Belgelendirme İşareti  
علامة المطابقة للمعايير الأوروبية CE

## Veraseal® 2 Full Face Non-Vented Mask

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The Veraseal 2 Full Face Non-Vented Mask is a disposable device intended to provide a patient interface for the application of noninvasive ventilation. The mask is to be used as an accessory to ventilators with exhalation valves, adequate alarms, and safety systems for ventilator failure, and which are intended to administer positive pressure ventilation. The mask is intended for SHORT TERM SINGLE PATIENT USE (maximum of 7 days) on adult patients (>30 kg) who are appropriate candidates for noninvasive ventilation in the hospital or institutional environment.

## NOTES

- ▶ The mask system does not contain PVC, DEHP or phthalates. If the patient has ANY reaction to any part of the mask system, discontinue use.
- ▶ This product is not manufactured with natural rubber latex.
- ▶ This mask requires a separate exhalation device.
- ▶ The mask does not include an anti-asphyxia valve.
- ▶ Verify that the mask is the right size.

## CAUTIONS

- ▶ In the U.S. Federal Law restricts this device for sale by or on the order of a physician.
- ▶ Do not expose the mask to temperatures above 140°F (60°C).
- ▶ The mask is intended for SINGLE PATIENT USE ONLY. It is intended for short term use (not more than 7 days) and must not be used on multiple persons. Do not sterilize or disinfect.
- ▶ Consult attending physician prior to using mask if the patient is using any medications or devices to remove deep secretions.

## WARNINGS

- ▶ This product should not be used if the patient is experiencing nausea, vomiting, taking a prescription drug that may cause vomiting, or if they are unable to remove the mask on their own.
- ▶ This mask is not suitable for providing life support ventilation.
- ▶ This mask must be used with a POSITIVE PRESSURE VENTILATOR WITH AN EXHALATION VALVE. The mask should not be used unless the positive pressure ventilator is turned on and operating properly. When the ventilator is functioning properly, the exhalation valve allows exhaled air to escape to ambient. When the ventilator is not operating, exhaled air may be re-breathed. Re-breathing of exhaled air can, in some circumstances, lead to suffocation.
- ▶ This mask must be used with VENTILATORS WITH EXHALATION VALVES, ADEQUATE ALARMS, AND SAFETY SYSTEMS FOR VENTILATOR FAILURE; use of this mask requires the same level of attention and assistance as the use of a tracheal tube.
- ▶ This mask does not include an anti-asphyxia valve that would allow a patient to breathe if the ventilator failed.
- ▶ If the patient has respiratory insufficiency, they should wear this mask ONLY when therapy is being delivered.
- ▶ Use this mask with baseline pressure of 3 cm H<sub>2</sub>O or higher.
- ▶ If oxygen is used with this mask, the oxygen flow must be turned off when the machine is not operating. Explanation of the warning: When the positive pressure device is not in operation, and the oxygen flow is left on, oxygen delivered into the device tubing may accumulate within the positive pressure device enclosure. Oxygen accumulated will create a risk of fire.
- ▶ At a fixed flow rate of supplemental oxygen, the inhaled oxygen concentration will vary, depending on pressure settings, patient breathing, and leak rate.
- ▶ No smoking or open flames such as candles may be used when oxygen is in use.
- ▶ This mask is not intended for patients without a spontaneous respiratory drive. This mask should not be used on patients who are uncooperative, obtunded, unresponsive, or unable to remove the mask. This mask may not be suitable for persons with the following conditions: impaired cardiac sphincter function, excessive reflux, impaired cough reflex, and hiatal hernia.
- ▶ Any unusual skin irritation, chest discomfort, shortness of breath, gastric distension, abdominal pain, belching or flatulence from ingested air, or severe headache experienced by the patient during or immediately after use should be immediately reported to the attending healthcare professional.
- ▶ Using a mask may cause tooth, gum or jaw soreness or aggravate an existing dental condition. If symptoms occur, consult a physician or dentist.
- ▶ Consult a physician if the patient experiences any of the following symptoms while using the mask: Drying of the eyes, eye pain, eye infections, or blurred vision. Consult an ophthalmologist if symptoms persist.
- ▶ Do not overtighten the headgear straps. Signs of overtightening include: excessive redness, sores, or bulging skin around the edges of the mask. Loosen headgear straps to alleviate symptoms.
- ▶ Significant unintentional leak may occur without properly fitting the mask.

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## GENERAL INFORMATION

The mask has a service life of 7 days. Please dispose of the mask and accessories in accordance with local and regional guidelines. This mask and all its accessories do not contain natural rubber latex. Per EU MDR 2017/745, please report serious incidents to manufacturer and EU Member State competent authority.

## CLEANING AND MAINTENANCE

This mask is disposable and is intended for single patient short term use (maximum of 7 days) only. The mask should be discarded if heavily soiled. The mask may be cleaned daily or as needed using a 70% v/v isopropyl alcohol swab or by wiping with a damp, soapy cloth. If using a damp cloth, use mild soap. Rinse thoroughly and allow to air dry, out of direct sunlight.

**NOTE:** The gel cushion, although silky to the touch, will not tear, leak, or detach from the mask under normal use. The gel cushion is not indestructible. It will tear, cut, and rip if subjected to mishandling or abuse.

## TECHNICAL INFORMATION

Resistance to air flow at 50 L/min: 0.10 cm H<sub>2</sub>O

Resistance to air flow at 100 L/min: 0.50 cm H<sub>2</sub>O

Dead Space (approx.): **Small** 155.3 ml **Medium** 181.7 ml **Large** 196.6 ml **Extra Large** 200.6 ml

Operating Temperature: Do not expose the mask to temperatures above 140°F (60°C).

## BEFORE USE

- Inspect mask daily or prior to each use. Replace the mask if parts are damaged, or if gel is exposed due to tears or punctures.

## MASK FITTING INSTRUCTIONS

- Detach one Velcro® strap by sliding the Velcro loop off the quick release hook, slip the headgear over the head and center the mask on the face (Figure 1).
  - Find the detached Velcro® strap and slide the Velcro loop back on the quick release hook (Figure 2).
  - Be sure the headgear straps are parallel and the headgear is fully extended and flat across the back of the head (Figure 3).
  - Adjust the bottom headgear straps by using the Velcro® tabs. After tightening, make sure the Active Headgear Connector™ can move freely. (Figure 4).
  - Adjust the upper headgear straps by using the Velcro® tabs (Figure 5).
- NOTE:** Do NOT overtighten the headgear. The Flexible Spacebar should not touch the forehead (Figure 5a).
- Adjust the bottom headgear straps to make sure the bottom of the mask is sealed comfortably against the patient's face.

**NOTE:** For optimal seal and comfort, mask should be snug, but not tight. Adjust the tension by detaching the Velcro straps, tightening or loosening the straps, and reattaching the Velcro.

- Connect the mask tubing assembly to the positive airway pressure device tubing and turn the device on to the prescribed pressure (Figure 6).

Velcro® is a registered trademark of Velcro Industries B.V.

Veraseal® is a trademark of Sleepnet Corporation.

Visit our website at [www.sleepnetmasks.com](http://www.sleepnetmasks.com).

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