V3 Full Face Vented Mask

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These instructions contain the following languages:

English, Spanish, German, French, Italian, Finnish, Turkish, Arabic, Bulgarian, Croatian, Czech, Danish, Dutch, Greek, Hungarian, Japanese, Norwegian, Polish, Portuguese, Romanian, Russian, Serbian, Slovak, Slovenian, Swedish and Hebrew

Sleepnet Corporation

5 Merrill Industrial Drive

Hampton, NH 03842 USA

Tel: (603) 758-6600

Fax: (603) 758-6699

Toll Free: 1-800-742-3646 (USA ONLY)

www.sleepnetmasks.com

90159 Rev D 2023-03


Figure 1


Figure 2


Figure 3


Figure 4


Figure 5


Figure 6


Figure 7a



Figure 7b

V3 Vented Mask

1. Forehead Pad
2. Adjustment Screw
3. Vent Holes
4. Airgel® Cushion
5. Elbow Assembly with AAV
6. Oxygen port
7. Headgear Connector
8. Quick Release Hook
9. Headgear

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
V3 Full Face Vented Mask

ENGLISH

CONTENTS: V3 Full Face Vented Mask with Headgear

Rx ONLY MADE IN U.S.A.

DOES NOT CONTAIN NATURAL RUBBER LATEX

INTENDED USE

The V3 Full Face Vented Mask is intended to be used with positive airway pressure devices, such as CPAP or bi-level, operating at or above 3 cm H₂O. The mask is to be used on adult patients (>30kg) for whom positive airway pressure therapy has been prescribed. The mask is intended for single-patient multi-use in the home, hospital or institutional environment.

NOTES

- The mask system does not contain PVC, DEHP, or phthalates.
- This product is not manufactured with natural rubber latex.
- If the patient has ANY reaction to any part of the mask system, discontinue use.
- An exhalation port is built into this mask so a separate exhalation port is not required.
- The mask includes an anti-asphyxia valve (air entrainment valve) to help reduce rebreathing of exhaled air in the event that the positive pressure device does not work properly.
- Verify that the mask is the right size by using the sizing guide.

CONTRAINDICATIONS

This product should not be used if the patient is experiencing nausea, vomiting, taking a prescription drug that may cause vomiting, or if they are unable to remove the mask on their own.

CAUTIONS

- U.S. Federal Law restricts this device for sale by or on the order of a physician.
- Do not expose the V3 Full Face Vented Mask to temperatures above 140°F (60°C).
- This mask is intended for single patient use only. It may be cleaned and used repeatedly on the same person but must not be used on multiple persons. Do not sterilize or disinfect.
- Consult a physician prior to using mask if using any medications or devices to remove deep secretions.

WARNINGS

- This mask is not suitable for providing life support ventilation.
- If the patient has respiratory insufficiency, they should wear this mask ONLY when therapy is being delivered.
- Do not use the V3 Full Face Vented Mask with baseline pressure less than 3 cm H₂O.
- If oxygen is used with this CPAP mask, the oxygen flow must be turned off when the CPAP machine is not operating. Explanation of the warning: When the CPAP device is not in operation, and the oxygen flow is left on, oxygen delivered into the ventilator tubing may accumulate within the CPAP machine enclosure. Oxygen accumulated in the CPAP machine enclosure will create a risk of fire.
- At a fixed rate of supplemental oxygen flow, the inhaled oxygen concentration will vary, depending on the pressure settings, patient breathing pattern, mask selection, and the leak rate.
- No smoking or open flames such as candles may be used when oxygen is in use.
- To minimize the risk of vomiting during sleep, avoid eating or drinking for three (3) hours before using the mask.
- Any unusual skin irritation, chest discomfort, shortness of breath, gastric distension, abdominal pain, belching or flatulence from ingested air, or severe headache experienced by the patient during or immediately after use should be immediately reported to a healthcare professional.
- Using a mask may cause tooth, gum or jaw soreness or aggravate an existing dental condition. If symptoms occur, consult a physician or dentist.
- Consult a physician if any of the following symptoms occur while using the mask: drying of the eyes, eye pain, eye infections, or blurred vision. Consult an ophthalmologist if symptoms persist.
- Do not overtighten the headgear straps. Signs of overtightening include: excessive redness, sores, or bulging skin around the edges of the mask. Loosen headgear straps to alleviate symptoms.
- Significant unintentional leak may occur without properly fitting the mask. Follow the MASK FITTING INSTRUCTIONS provided.
- The elbow and anti-asphyxia valve assembly have specific safety functions. The mask should not be used if the valve flapper is damaged, distorted, or torn.
- At low CPAP pressures, the flow through the exhalation ports may be inadequate to clear all exhaled gas from the tubing. Some rebreathing may occur.
- This mask should be used with the positive airway pressure device recommended by a physician or respiratory therapist. The mask should not be used unless the positive airway pressure device is turned on and operating properly. Do not block or try to seal the exhalation port (vent holes) or anti-asphyxia valve (air entrainment valve). The vent holes allow a continuous flow of air out of the mask. When the positive airway pressure device is turned on and functioning properly, new air from the device flushes the exhaled air out through the mask's vent holes. When the positive airway pressure device is not operating, exhaled air may be rebreathed. Rebreathing of exhaled air can, in some circumstances, lead to suffocation. This warning applies to most models of CPAP machines and masks.
- This mask is not intended for patients without a spontaneous respiratory drive. This mask should not be used on patients who are uncooperative, obtunded, unresponsive, or unable to remove the mask. This mask may not be suitable for persons with the following conditions: impaired cardiac sphincter function, excessive reflux, impaired cough reflex, and hiatal hernia.

GENERAL INFORMATION

This mask has a non-transferable three month limited warranty for manufacturer's defects from the date of purchase by the initial consumer. If a mask used under normal conditions should fail, Sleepnet will replace the mask. For additional information about the Sleepnet Warranty, please visit <http://www.sleepnetmasks.com>. Per EU MDR 2017/745, please report serious incidents to manufacturer and EU Member State competent authority.

NOTE: The mask has a service life of 6 months.

CLEANING AND MAINTENANCE

Cleaning is recommended on a daily basis. Wash the mask and components in warm water using a mild soap such as Ivory®. Rinse the mask thoroughly after washing. Visually verify that the vent holes are clear and allow the mask to air dry. Washing your face prior to fitting the mask will remove excess facial oils and help extend the life of the mask.

NOTE: Inspect mask and valve flapper prior to each use. Replace the mask if parts are damaged or gel is exposed due to tear or punctures.

NOTE: The gel cushion, although silky to touch, will not tear, leak or detach from the mask under normal use and following recommended cleaning and maintenance protocols. However, the gel cushion is not indestructible. It will tear, cut and rip if subjected to mishandling or abuse. Do not clean in a dishwashing machine. The warranty does not cover these types of product failures. Therefore, please handle your new V3 Full Face Vented Mask with care.

TECHNICAL INFORMATION

Resistance to air flow at 50 L/min: 0.23 cm H₂O; at 100 L/min: 0.48 cm H₂O

DEADSPACE VOLUME (mL approx.) **Small** 184 mL **Medium** 193 mL **Large** 209 mL **Extra Large** 220 mL

Resistance of the anti-asphyxia valve in Single Fault Condition: Inspiratory: 0.8 H₂O per L/s. Expiratory: 0.8 H₂O per L/s

A-weighted Sound Pressure Level at 1 m distance: 29.98 (A); A-weighted Sound Power Level: 32.99 (A)

AAV Open to Atmosphere Pressure: 1.45 cm H₂O, AAV Close to Atmosphere Pressure: 1.80 cm H₂O

Do not expose the mask to temperatures above 140°F (60°C).

BEFORE USE

- Inspect mask daily or prior to each use. Replace the mask if vent holes are blocked, if parts are damaged, or if gel is exposed due to tears or punctures.
- Check the anti-asphyxia valve. With the positive pressure device turned off, verify that the valve flapper is positioned so that room air can flow through the large opening in the valve (Figure 7a). Turn on the system; the flapper should close, and system air should flow through the mask (Figure 7b). If the flapper does not close or does not function properly, replace the mask. Do not block the opening on the anti-asphyxia valve or the vent holes. Ensure that the valve is not blocked with secretions and that the flapper is dry.

V3 VENTED MASK FITTING INSTRUCTIONS

1. Verify the mask is the right size. When using the sizing guide ensure mouth is slightly open. Line the top of the guide to the bridge of the nose and select the smallest size that does not restrict the nose or mouth. The appropriate size may be smaller than expected due to the design.
2. Lay headgear flat with light gray side facing up with the shorter straps at the top. Place mask in the center of headgear and thread 3 straps as shown (Figure 1).
3. Seat the bottom cushion of the mask at the crease of the chin with your mouth slightly open and tip the mask up to contact the nose (Figure 2). Turn the forehead adjustment screw so the forehead pad lightly touches your forehead (Figure 6). The mask may sit lower on the nose than expected, depending on facial features.
4. Slip the headgear over the head. Connect the detached strap over the opening of the Quick Release slot (Figure 3). The headgear should be flat across the back of the head.
5. Adjust the bottom headgear straps by using the hook and loop fastener tabs (Figure 4).
6. Adjust the upper headgear straps: detach the hook and loop fastener tabs and pull straps forward to center mask over face. Ensure headgear is centered by aligning the two straps. Then pull the straps back, adjust strap tension and affix tabs. (Figure 5)
7. Turn the forehead pad adjustment screw to adjust the tension until comfortable. If necessary, readjust top straps for a comfortable fit.
8. If the mask is uncomfortable repeat steps 5-8.

NOTE:

- I. If breathing through your nose or mouth is still constricted, contact your healthcare professional.
- II. The mask should rest comfortably against the face. If there is excessive pressure on the bridge of the nose, loosen the top straps, slowly turn the adjustment screw counter-clockwise to relieve the pressure on the bridge (Figure 6). Then readjust the top straps and affix the tabs.
9. Connect the mask to the positive pressure device tubing and breathe normally.
10. If leaks are noticed around the chin or the cheeks, slightly tighten the bottom straps.
11. If leaks are noticed at the bridge of the nose, slowly turn the adjustment screw clockwise to slightly move the pad away from the forehead (Figure 6). Tighten the top straps slightly until leaks are eliminated.

NOTE:

- I. DO NOT overtighten any straps. Overtightening can worsen leaks. For optimal seal and comfort, mask should be snug, but not tight.

QUICK RELEASE OF MASK AND HEADGEAR

- Slide lower right headgear strap off the quick release slot and remove mask.

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Medical Product Services GmbH