Explanation of Symbols Which Appear on the Device Labeling
Explicación de los símbolos que aparecen en la etiqueta del dispositivo
Erklärung der Symbole auf dem Etikett dieses Gerätes
Explication des symboles qui apparaissent sur l'étiquette de l'appareil
Spiegazione dei simboli che figurano sulle etichette del dispositivo
Laitteen etiketissä olevien symbolien selitykset
装置标签上的符号说明
装置のラベルに表示される記号の説明

Caution or Warning
Precaución o advertencia
Vorsichts- oder Warnhinweise
Avertissement ou mise en garde
Attenzione o avvertenza
Huomio tai varoitus
注意事項或警告
注意・警告

Refer to Instructions for Use
Consulte las instrucciones de uso.
Siehe Gebrauchsanleitung
Se reporter au Mode d'emploi
Consultare le Istruzioni per l'uso
Katso käyttöohjeet
请参阅使用说明
使用説明書を参照

Manufacturer
Fabricante
Hersteller
Fabricant
Fabbicante
Valmistaja
制造商
製造元

CE Certification Mark
Marca de certificación CE
CE Zertifizierungsmarke
Marquage CE
Marcatura CE
CE-sertifiointimerkki
CE 认证标记
CE認証マーク

Part number
Referencia
Teilenummer
Référéce
Codice articolo
Osanumero
部件号
部品番号

Lot number
Número de lote
Losnummer
Numéro du lot
Numero di lotto
Eränumero
批号
ロット番号

This Product is not manufactured with Natural Ruber Latex.
Este producto no se fabrica con látex de caucho natural.
Dieses Produkt ist nicht mit Naturkautschuklatex hergestellt
Ce produit est fabriqué sans latex de caoutchouc naturel.
Il prodotto non contiene lattice di gomma naturale.
Tuotteen valmistuksessa ei ole käytetty luonnonkumilateksia.
本产品不含天然胶乳。
この製品には天然ゴムラテックスは使用されていません。

Date of Manufacture
Fecha de fabricación
Herstellungsdatum
Date de fabrication
Data di produzione
Valmistuspäivämäärä
制造日期
製造年月日

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This page has been left blank intentionally

Veraseal 2 Full Face Non-Vented AAV Mask

Contents: Veraseal 2 Non-Vented AAV Mask with headgear

Made in the U.S.A. of U.S. and imported parts

Rx only

INTENDED USE

The Veraseal 2 Full Face Non-Vented AAV Mask is a disposable device intended to provide a patient interface for the application of noninvasive ventilation. The mask is to be used as an accessory to ventilators with exhalation valves, adequate alarms, and safety systems for ventilator failure, and which are intended to administer positive pressure ventilation. The mask is intended for SHORT TERM SINGLE PATIENT USE (maximum of 7 days) on adult patients (>30 kg) who are appropriate candidates for noninvasive ventilation in the hospital or institutional environment.

NOTES

- ▶ The mask system does not contain PVC, DEHP, or phthalates. If the patient has ANY reaction to any part of the mask system, discontinue use.
- ▶ This product is not manufactured with natural rubber latex.
- ▶ This mask requires a separate exhalation device.
- ▶ The mask includes an anti-asphyxia valve to help reduce rebreathing of exhaled air in the event that the positive pressure device does not work properly.
- ▶ Verify that the mask is the right size.

CAUTIONS

- ▶ In the U.S. Federal Law restricts this device for sale by or on the order of a physician.
- ▶ Do not expose the mask to temperatures above 140°F (60°C).
- ▶ The mask is intended for SINGLE PATIENT USE ONLY. It is intended for short term use (not more than 7 days) and must not be used on multiple persons. Do not sterilize or disinfect.
- ▶ Consult attending physician prior to using mask if the patient is using any medications or devices to remove deep secretions.

WARNINGS

- ▶ This product should not be used if the patient is experiencing nausea, vomiting, taking a prescription drug that may cause vomiting, or if they are unable to remove the mask on their own.
- ▶ This mask is not suitable for providing life support ventilation.
- ▶ The elbow and anti-asphyxia valve assembly have specific safety functions. The mask should not be used if the valve is damaged. The mask should be discarded if the valve is damaged, distorted, or torn.
- ▶ If the patient has respiratory insufficiency, they should wear this mask ONLY when therapy is being delivered.
- ▶ Use this mask with baseline pressure of 3 cm H₂O or higher.
- ▶ This mask must be used with a POSITIVE PRESSURE VENTILATOR WITH AN EXHALATION VALVE. The mask should not be used unless the positive pressure ventilator is turned on and operating properly. When the ventilator is functioning properly, the exhalation valve allows exhaled air to escape to ambient. When the ventilator is not operating, exhaled air may be re-breathed. Re-breathing of exhaled air can, in some circumstances, lead to suffocation.
- ▶ This mask must be used with VENTILATORS WITH EXHALATION VALVES, ADEQUATE ALARMS, AND SAFETY SYSTEMS FOR VENTILATOR FAILURE; use of this mask requires the same level of attention and assistance as the use of a tracheal tube.
- ▶ If oxygen is used with this CPAP mask, the oxygen flow must be turned off when the CPAP machine is not operating. Explanation of the warning: When the CPAP device is not in operation, and the oxygen flow is left on, oxygen delivered into the device tubing may accumulate within the CPAP machine enclosure. Oxygen accumulated in the CPAP machine enclosure will create a risk of fire.
- ▶ At a fixed flow rate of supplemental oxygen, the inhaled oxygen concentration will vary, depending on pressure settings, patient breathing, and leak rate.
- ▶ No smoking or open flames such as candles may be used when oxygen is in use.
- ▶ This mask is not intended for patients without a spontaneous respiratory drive. This mask should not be used on patients who are uncooperative, obtunded, unresponsive, or unable to remove the mask. This mask may not be suitable for persons with the following conditions: impaired cardiac sphincter function, excessive reflux, impaired cough reflex, and hiatal hernia.
- ▶ Any unusual skin irritation, chest discomfort, shortness of breath, gastric distension, abdominal pain, belching or flatulence from ingested air, or severe headache experienced by the patient during or immediately after use should be immediately reported to the attending healthcare professional.

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- ▶ Using a mask may cause tooth, gum or jaw soreness or aggravate an existing dental condition. If symptoms occur, consult a physician or dentist.
- ▶ Consult a physician if the patient experiences any of the following symptoms while using the mask: Drying of the eyes, eye pain, eye infections, or blurred vision. Consult an ophthalmologist if symptoms persist.
- ▶ Do not overtighten the headgear straps. Signs of overtightening include: excessive redness, sores, or bulging skin around the edges of the mask. Loosen headgear straps to alleviate symptoms.
- ▶ Significant unintentional leak may occur without properly fitting the mask.

GENERAL INFORMATION

The mask has a service life of 7 days. Please dispose of the mask and accessories in accordance with local and regional guidelines. This mask and all its accessories do not contain natural rubber latex.

CLEANING AND MAINTENANCE

This mask is disposable and is intended for single patient short term use (maximum of 7 days) only. The mask should be discarded if heavily soiled. The mask may be cleaned daily or as needed using a 70% v/v isopropyl alcohol swab or by wiping with a damp, soapy cloth. If using a damp cloth, use mild soap. Rinse thoroughly and allow to air dry, out of direct sunlight.

NOTE: The gel cushion, although silky to the touch, will not tear, leak, or detach from the mask under normal use. The gel cushion is not indestructible. It will tear, cut, and rip if subjected to mishandling or abuse.

TECHNICAL INFORMATION

Resistance to air flow at 50 L/min: 0.07 cm H₂O
Resistance to air flow at 100 L/min: 0.39 cm H₂O
Dead Space (approx.): **Small** 155.3 ml **Medium** 181.7 ml **Large** 196.6 ml **Extra Large** 200.6 ml
Anti-asphyxia valve: Open-to-atmosphere pressure: 2.1 cm H₂O
Closed-to-atmosphere pressure: 2.28 cm H₂O
Operating Temperature: Do not expose the mask to temperatures above 140°F (60°C).

BEFORE USE

- ▶ Inspect mask daily or prior to each use. Replace the mask if parts are damaged, or if gel is exposed due to tears or punctures.
- ▶ Check the anti-asphyxia valve. With the positive pressure device turned off, verify that the valve flapper is positioned so that room air can flow through the large opening in the valve (Figure 1a). Turn on the system; the flapper should close, and system air should flow through the mask (Figure 1b). If the flapper does not close or does not function properly, replace the mask. Do not block the opening on the anti-asphyxia valve or the vent holes. Ensure that the valve is not blocked with secretions and that the flapper is dry.

MASK FITTING INSTRUCTIONS

1. Detach one Velcro® strap by sliding the Velcro loop off the quick release hook, slip the headgear over the head and center the mask on the face (Figure 2).
2. Find the detached Velcro® strap and slide the Velcro loop back on the quick release hook (Figure 3).
3. Be sure the headgear straps are parallel and the headgear is fully extended and flat across the back of the head (Figure 4).
4. Adjust the bottom headgear straps by using the Velcro® tabs. After tightening, make sure the Active Headgear Connector™ can move freely. (Figure 5).
5. Adjust the upper headgear straps by using the Velcro® tabs (Figure 6).

NOTE: Do NOT overtighten the headgear. The Flexible Spacebar should not touch the forehead (Figure 6a).

6. Adjust the bottom headgear straps to make sure the bottom of the mask is sealed comfortably against the patient's face.

NOTE: For optimal seal and comfort, mask should be snug, but not tight. Adjust the tension by detaching the Velcro straps, tightening or loosening the straps, and reattaching the Velcro.

7. Connect the mask tubing assembly to the positive airway pressure device tubing and turn the device on to the prescribed pressure (Figure 7).

Velcro® is a registered trademark of Velcro Industries B.V.

Veraseal® is a trademark of Sleepnet Corporation.

Visit our website at www.sleepnetmasks.com.

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