Manufacturer
Fabricante
Hersteller
Fabricant
Fabbriçante
Valmistaja
製造商
製造元
Üretici
Изготовитель
Producător
Tillverkare
Производитель
Proizvodjač
Výrobce
Producent
Fabrikant
Κατασκευαστής
Gyártó
Produsent
Producent
Fabricante
Proizvodjač
Výrobca
Izdelovalec
الجهة المصنعة

Date of Manufacture
Fecha de fabricación
Herstellungsdatum
Date de fabrication
Data di produzione
Valmistuspäivämäärä
製造日期
製造年月日
Üretim Tarihi
Дата изготовления
Data fabricației
Tillverkningsdatum
Дата на производство
Datum proizvodnje
Datum výroby
Fremstillingsdato
Productiedatum
Ημερομηνία κατασκευής
A gyártás dátuma
Produksjonsdato
Data produkci
Data de fabrico
Datum proizvodnje
Datum výroby
Datum izdelave
تاريخ التصنيع

CE Certification Mark
Marca de certificación CE
CE Zertifizierungsmarke
Marquage CE
Marcatura CE
CE-sertifiointimerkki
CE 认证标记
CE 認證マーク
CE Belgeleendirme İşareti
Знак сертификации CE
Marcă de certificare CE
CE-märkning
CE сертификационная маркировка
Certifikacijska oznaka CE
Označení CE certifikace
CE-mærke
CE-markering
Σήμα πιστοποίησης CE
CE tanúsítvány-jelölés
CE-sertifiseringsmerke
Znak CE
Marca de certificação CE
CE oznaka
Certifikačná známka CE
Oznaka CE
CE علامة المطابقة للمعايير الأوروبية

Part number
Referencia
Teilenummer
Référence
Codice articolo
Osnumero
部件号
部品番号
Parti numarası
Номер по каталогу
Număr piesă
Art.nr
Номер на част
Broj dijela
Číslo dílu
Delnummer
Onderdeelnnummer
Αριθμός εξαρτήματος
Cikkszám
Delenummer
Numer katalogowy
Número de peça
Serijski broj
Číslo dielu
Številka dela
رقم القطعة


Medical Device
Producto sanitario
Medizinisches Gerät
Dispositif médical
Dispositivo medico
Lääkinnällinen laite
医疗器械
医療機器
Medikal Cihaz
Медицинское устройство
Medicinski pripomoček
Medicinsk utrustning
Медицинско изделие
Medicinski proizvod
Zdravotnický prostriedek
Medicinsk enhed
Medisch hulpmiddel
Ιατρική Συσκευή
Orvostechnikai eszköz
Medisinsk enhet
Urządzenie medyczne
Dispositivo médico
Medicinsko sredstvo
Zdravotnícka pomôcka
Medicinski pripomoček
جهاز طبي

Use by Date
Fecha de caducidad
Verfallsdatum
Date limite d'utilisation
Da utilizzare entro la data indicata
Viimeinen käyttöpäivä
此日期前使用
使用期限
Son Kullanma Tarihi
Использовать до указанной даты
Datum upotrebe
Bör användas före
Крайна дата на използване
Upotrijebiti do datuma
Datum použitelnosti
Sidste anvendelsesdato
Uiterste houdbaarheidsdatum
Ημερομηνία λήξης
Felhasználhatósági idő
Brukes innen
Data przydatności do użycia
Prazo de validade
Rok upotrebe
Datum trvanlivosti
Rok uporabnosti
يستخدم قبل التاريخ

Importer/ Importador / Importeur/ Importateur/ Importatore/ Maahantuoja/ 进口商 / 輸入業者 / İthalatçı/
Импортер/ Importator/ Importör/ Вноситель/ Uvoznik/ Importér/ Importor/ Importeur/ Εισαγωγέας/
Importör/ Importør/ Importer/ Importador/ Uvoznik/ Dovožca/ Uvoznik/ المستورد

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
VERASEAL® 2 FULL FACE VENTED MASK

Contents: Veraseal 2 Vented Mask with headgear

Made in the U.S.A. of U.S. and imported parts

Rx only

INTENDED USE

The Veraseal® 2 Full Face Vented Mask is a disposable mask intended to be used with positive airway pressure devices, such as CPAP or bi-level, operating at or above 3 cm H₂O. The mask is to be used on adult patients (>30kg) for whom positive airway pressure therapy has been prescribed. The mask is intended for SHORT TERM SINGLE PATIENT USE (maximum of 7 days) in the hospital or institutional environment.

NOTES

- The mask system does not contain PVC, DEHP, or phthalates. If the patient has ANY reaction to any part of the mask system, discontinue use.
- This product is not manufactured with natural rubber latex.
- An exhalation port is built into this mask so a separate exhalation port is not required.
- The mask includes an anti-asphyxia valve to help reduce rebreathing of exhaled air in the event that the positive pressure device does not work properly.
- Verify that the mask is the right size.

CONTRAINDICATION

This product should not be used if the patient is experiencing nausea, vomiting, taking a prescription drug that may cause vomiting, or if they are unable to remove the mask on their own.

CAUTIONS

- In the U.S. Federal Law restricts this device for sale by or on the order of a physician.
- Do not expose the mask to temperatures above 140°F (60°C).
- The mask is intended for SINGLE PATIENT USE ONLY. It is intended for short term use (not more than 7 days) and must not be used on multiple persons. Do not sterilize or disinfect.
- Consult attending physician prior to using mask if the patient is using any medications or devices to remove deep secretions.

WARNINGS

- This mask is not suitable for providing life support ventilation.
- The elbow and anti-asphyxia valve assembly have specific safety functions. The mask should not be used if the valve is damaged. The mask should be discarded if the valve is damaged, distorted, or torn.
- If the patient has respiratory insufficiency, they should wear this mask ONLY when therapy is being delivered.
- At low pressures, the flow through the exhalation ports may be inadequate to clear all exhaled gas from the tubing. Some rebreathing may occur.
- Use this mask with baseline pressure of 3 cm H₂O or higher.
- Do not block or try to seal the vent holes or anti-asphyxia valve. The vent holes allow a continuous flow of air out of the mask.
- The mask should not be used unless the positive airway pressure device is turned on and operating properly. Explanation of the warning: When the positive airway pressure device is turned on and functioning properly, new air from the device flushes the exhaled air out through the mask's vent holes. When the positive airway pressure device is not operating, exhaled air may be rebreathed. Rebreathing of exhaled air can, in some circumstances, lead to suffocation. This warning applies to most models of CPAP machines and masks.
- If oxygen is used with this CPAP mask, the oxygen flow must be turned off when the CPAP machine is not operating. Explanation of the warning: When the CPAP device is not in operation, and the oxygen flow is left on, oxygen delivered into the device tubing may accumulate within the CPAP machine enclosure. Oxygen accumulated in the CPAP machine enclosure will create a risk of fire.
- At a fixed flow rate of supplemental oxygen, the inhaled oxygen concentration will vary, depending on pressure settings, patient breathing, and leak rate.
- No smoking or open flames such as candles may be used when oxygen is in use.
- This mask is not intended for patients without a spontaneous respiratory drive. This mask should not be used on patients who are uncooperative, obtunded, unresponsive, or unable to remove the mask. This mask may not be suitable for persons with the following conditions: impaired cardiac sphincter function, excessive reflux, impaired cough reflex, and hiatal hernia.
- Any unusual skin irritation, chest discomfort, shortness of breath, gastric distension, abdominal pain, belching or flatulence from ingested air, or severe headache experienced by the patient during or immediately after use should be immediately reported to the attending healthcare professional.
- Using a mask may cause tooth, gum or jaw soreness or aggravate an existing dental condition. If symptoms occur, consult a physician or dentist.
- Consult a physician if the patient experiences any of the following symptoms while using the mask: Drying of the eyes, eye pain, eye infections, or blurred vision. Consult an ophthalmologist if symptoms persist.
- Do not overtighten the headgear straps. Signs of overtightening include: excessive redness, sores, or bulging skin around the edges of the mask. Loosen headgear straps to alleviate symptoms.
- Significant unintentional leak may occur without properly fitting the mask.

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GENERAL INFORMATION

The mask has a service life of 7 days. Please dispose of the mask and accessories in accordance with local and regional guidelines. This mask and all its accessories do not contain natural rubber latex. Per EU MDR 2017/745, please report serious incidents to manufacturer and EU Member State competent authority.

CLEANING AND MAINTENANCE

This mask is disposable and is intended for single patient short term use (maximum of 7 days) only. The mask should be discarded if heavily soiled. The mask may be cleaned daily or as needed using a 70% v/v isopropyl alcohol swab or by wiping with a damp, soapy cloth. If using a damp cloth, use mild soap. Rinse thoroughly and allow to air dry, out of direct sunlight.

NOTE: The gel cushion, although silky to the touch, will not tear, leak, or detach from the mask under normal use. The gel cushion is not indestructible. It will tear, cut, and rip if subjected to mishandling or abuse.

TECHNICAL INFORMATION

Resistance to air flow at 50 L/min: 0.07 cm H₂O

Resistance to air flow at 100 L/min: 0.39 cm H₂O

DEAD SPACE (approx.): **Small** 155.3 ml **Medium** 181.7 ml **Large** 196.6 ml **Extra Large** 200.6 ml

Anti-asphyxia valve: Open-to-atmosphere pressure: 2.1 cm H₂O

Closed-to-atmosphere pressure: 2.28 cm H₂O

Operating Temperature: Do not expose the mask to temperatures above 140°F (60°C).

BEFORE USE

- Inspect mask daily or prior to each use. Replace the mask if vent holes are blocked, if parts are damaged, or if gel is exposed due to tears or punctures.
- Check the anti-asphyxia valve. With the positive pressure device turned off, verify that the valve flapper is positioned so that room air can flow through the large opening in the valve (Figure 1a). Turn on the system; the flapper should close, and system air should flow through the mask (Figure 1b). If the flapper does not close or does not function properly, replace the mask. Do not block the opening on the anti-asphyxia valve or the vent holes. Ensure that the valve is not blocked with secretions and that the flapper is dry.

MASK FITTING INSTRUCTIONS

1. Detach one Velcro® strap by sliding the Velcro loop off the quick release hook, slip the headgear over the head and center the mask on the face (Figure 2).
2. Find the detached Velcro® strap and slide the Velcro loop back on the quick release hook (Figure 3).
3. Be sure the headgear straps are parallel and the headgear is fully extended and flat across the back of the head (Figure 4).
4. Adjust the bottom headgear straps by using the Velcro® tabs. After tightening, make sure the Active Headgear Connector™ can move freely. (Figure 5).
5. Adjust the upper headgear straps by using the Velcro® tabs (Figure 6).

NOTE: Do NOT overtighten the headgear. The Flexible Spacebar should not touch the forehead (Figure 6a).

6. Adjust the bottom headgear straps to make sure the bottom of the mask is sealed comfortably against the patient's face.

NOTE: For optimal seal and comfort, mask should be snug, but not tight. Adjust the tension by detaching the Velcro straps, tightening or loosening the straps, and reattaching the Velcro.

7. Connect the mask tubing assembly to the positive airway pressure device tubing and turn the device on to the prescribed pressure (Figure 7).

Velcro® is a registered trademark of Velcro Industries B.V.

Veraseal® is a trademark of Sleepnet Corporation.

Visit our website at www.sleepnetmasks.com.

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