Mojo® 2 Full Face Vented Mask

- Mascarilla facial completa con ventilación Mojo® 2
Mojo® 2 Ganzgesichtsmaske (belüftet)
Masque ventilé Mojo® couvrant à la fois le nez et la bouche (2 unités)
Maschera ventilata integrale Mojo® 2
Mojo® 2 Geventileerd volgezichtsmasker
Mojo® 2 Αεριζόμενη Στοματινική Μάσκα
Mojo® 2-ansigtsiltmaske med ventilation
Mojo® 2 Ventilrad helansiktsmask
Máscara integral ventilada Mojo® 2
Pełnotwarzowa maska z otworami wentylacyjnymi Mojo® 2
Mojo® 2 aukollinen kokomaski
Mojo® 2 Havalandırılmalı Tam Yüz Maske
كمامة مُهَوّاة لكامل الوجه Mojo® 2
Mojo® 2 маска за цяло лице с вентилация
Ventilirana maska za cijelo lice Mojo® 2
Celoobličejová maska Mojo® 2 s větracími otvory
Mojo® 2 szellőző teljes arcmaszk
Mojo® 2 通气孔付きフルフェイスマスク
Mojo® 2 ventilert helansiktsmaske
Mască facială integrală Mojo® 2 cu ventilație
Вентилируемая полнолицевая маска Mojo® 2
Mojo® 2 maska za celo lice sa ventilacijom
Celotvárová ventilovaná maska Mojo® 2
Ventilirana celoobrazna maska Mojo® 2
חסימת פנים מלאה בעלת חורי אוורור - Mojo® 2

These instructions contain the following languages:

English, Spanish, German, French, Italian, Dutch, Greek, Danish, Swedish, Portuguese,
Polish, Finnish, Turkish, Arabic, Bulgarian, Croatian, Czech, Hungarian, Japanese, Norwegian,
Romanian, Russian, Serbian, Slovak, Slovenian and Hebrew

Sleepnet Corporation

5 Merrill Industrial Drive

Hampton, NH 03842 USA

Tel: (603) 758-6600

Fax: (603) 758-6699

Toll Free: 1-800-742-3646 (USA ONLY)

www.sleepnetmasks.com


Figure 1


Figure 2


Figure 3


Figure 4


Figure 5


Figure 6


Figure 7a



Figure 7b

Mojo® 2 Vented Mask

1. Forehead Pad
2. Adjustment Screw
3. Vent Holes
4. Airgel® Cushion
5. Elbow Assembly with AAV
6. Oxygen port
7. Magnetic clips
8. Pliable ring
9. Headgear

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
CONTENTS: Mojo® 2 Full Face Vented Mask with Headgear
Assembled in USA with US and imported parts
Rx Only

INTENDED USE

The Mojo® 2 Full Face Vented Mask is intended to be used with positive airway pressure devices, such as CPAP or bi-level, operating at or above 3 cm H₂O. The mask is to be used on adult patients (>30kg) for whom positive airway pressure therapy has been prescribed. The mask is intended for single-patient use in the home, hospital or institutional environment.

NOTES

- The mask system does not contain DEHP or phthalates.
- This product is not manufactured with natural rubber latex.
- If the patient has ANY reaction to any part of the mask system, discontinue use.
- An exhalation port is built into this mask so a separate exhalation port is not required.
- The mask includes an anti-asphyxia valve (air entrainment valve) to help reduce rebreathing of exhaled air in the event that the positive pressure device does not work properly.
- Verify that the mask is the right size by using the sizing guide.

CONTRAINDICATIONS

- This product should not be used if the patient is experiencing nausea, vomiting, taking a prescription drug that may cause vomiting, or if they are unable to remove the mask on their own.
- Do not use this mask if you or anyone (example: household members, bed partners, caregivers, etc.) in close physical contact with your mask has an active medical implant or metallic implant that will interact with magnets. Implant examples include, but are not limited to, pacemakers, implantable cardioverter defibrillators (ICD), neurostimulators, aneurysm clips, metallic stents, ocular implants, insulin/infusion pumps, cerebral spinal fluid (CSF) shunts, embolic coils, metallic splinter, implants to restore hearing or balance with implanted magnets (such as cochlear implants), flow disruption devices, contact lenses with metal, dental implants, metallic cranial plates, screws, burr hole covers, bone substitute device, magnetic metallic implants/electrodes/valves placed in upper limbs, torso, or higher, etc. If you have any questions regarding the implant, consult your physician or the manufacturer of your implant.

CAUTIONS

- U.S. Federal Law restricts this device for sale by or on the order of a physician.
- Do not expose the Mojo® 2 Full Face Vented Mask to temperatures above 122°F (50°C).
- This mask is intended for single patient use only. It may be cleaned and used repeatedly on the same person, but must not be used on multiple persons. Do not sterilize or disinfect.
- Consult a physician prior to using mask if using any medications or devices to remove deep secretions.

WARNINGS

- This mask is not suitable for providing life support ventilation
- If the patient has respiratory insufficiency, they should wear this mask ONLY when therapy is being delivered.
- Do not use the Mojo® 2 Full Face Vented Mask with baseline pressure less than 3 cm H₂O.
- If oxygen is used with this CPAP mask, the oxygen flow must be turned off when the CPAP machine is not operating. Explanation of the warning: When the CPAP device is not in operation, and the oxygen flow is left on, oxygen delivered into the ventilator tubing may accumulate within the CPAP machine enclosure. Oxygen accumulated in the CPAP machine enclosure will create a risk of fire.
- At a fixed rate of supplemental oxygen flow, the inhaled oxygen concentration will vary, depending on the pressure settings, patient breathing pattern, mask selection, and the leak rate.
- No smoking or open flames such as candles may be used when oxygen is in use.
- To minimize the risk of vomiting during sleep, avoid eating or drinking for three (3) hours before using the mask.
- Any unusual skin irritation, chest discomfort, shortness of breath, gastric distension, abdominal pain, belching or flatulence from ingested air, or severe headache experienced by the patient during or immediately after use should be immediately reported to a healthcare professional.
- Using a mask may cause tooth, gum or jaw soreness or aggravate an existing dental condition. If symptoms occur, consult a physician or dentist.
- Consult a physician if any of the following symptoms occur while using the mask: Drying of the eyes, eye pain, eye infections, or blurred vision. Consult an ophthalmologist if symptoms persist.
- Do not overtighten the headgear straps. Signs of overtightening include: excessive redness, sores, or bulging skin around the edges of the mask. Loosen headgear straps to alleviate symptoms.
- Significant unintentional leak may occur without properly fitting the mask. Follow the MASK FITTING INSTRUCTIONS provided.
- The elbow and anti-asphyxia valve assembly have specific safety functions. The mask should not be used if the valve flapper is damaged, distorted, or torn.
- At low CPAP pressures, the flow through the exhalation ports may be inadequate to clear all exhaled gas from the tubing. Some rebreathing may occur.
- This mask should be used with the positive airway pressure device recommended by a physician or respiratory therapist. The mask should not be used unless the positive airway pressure device is turned on and operating properly. Do not block or try to seal the exhalation port (vent holes) or anti-asphyxia valve (air entrainment valve). The vent holes allow a continuous flow of air out of the mask. When the positive airway pressure device is turned on and functioning properly, new air from the device flushes the exhaled air out through the mask's vent holes. When the positive airway pressure device is not operating, exhaled air may be rebreathed. Rebreathing of exhaled air can, in some circumstances, lead to suffocation. This warning applies to most models of CPAP machines and masks.
- This mask is not intended for patients without a spontaneous respiratory drive. This mask should not be used on patients who are uncooperative, obtunded, unresponsive, or unable to remove the mask. This mask may not be suitable for persons with the following conditions: impaired cardiac sphincter function, excessive reflux, impaired cough reflex, and hiatal hernia.

- Magnets are used in the mask and headgear clips with a field strength of 380mT. With the exception of the devices identified in the contraindication, ensure that the mask is kept at least 6 inches (approx.16 cm) away from any other medical implants or medical devices that can be impacted by the magnetic fields to avoid possible effects from localized magnetic fields. This applies to you or anyone in close physical contact with your mask.

GENERAL INFORMATION

This mask has a non-transferable three month limited warranty for manufacturer's defects from the date of purchase by the initial consumer. If a mask used under normal conditions should fail, Sleepnet will replace the mask. For additional information about the Sleepnet Warranty, please visit <http://www.sleepnetmasks.com>. Per EU MDR 2017/745, please report serious incidents to manufacturer and EU Member State competent authority.

NOTE: The mask has a service life of 6 months.

CLEANING AND MAINTENANCE

Cleaning is recommended on a daily basis. Wash the mask and components in warm water using a mild soap such as Ivory®. Rinse the mask thoroughly after washing. Visually verify that the vent holes are clear and allow the mask to air dry. Washing your face prior to fitting the mask will remove excess facial oils and help extend the life of the mask.

NOTE: Inspect mask and valve flapper prior to each use. Replace the mask if parts are damaged or gel is exposed due to tear or punctures.

NOTE: The gel cushion, although silky to touch, will not tear, leak or detach from the mask under normal use and following recommended cleaning and maintenance protocols. However, the gel cushion is not indestructible. It will tear, cut and rip if subjected to mishandling or abuse. Do not clean in a dishwashing machine. The warranty does not cover these types of product failures. Therefore, please handle your new Mojo 2 Full Face Vented Mask with care.

TECHNICAL INFORMATION

Resistance to air flow at 50 L/min: 0.25 cm H₂O; at 100 L/min: 0.50 cm H₂O

DEADSPACE VOLUME (mL approx.) **Small** 179 mL **Medium** 190 mL **Large** 207 mL **Extra Large** 220 mL

Resistance of the anti-asphyxia valve in Single Fault Condition: Inspiratory: 0.8 H₂O per L/s. Expiratory: 0.8 H₂O per L/s

A-weighted Sound Pressure Level at 1 m distance: 30.07 (A); A-weighted Sound Power Level: 33.09 (A)

AAV Open to Atmosphere Pressure: 1.35 cm H₂O, AAV Close to Atmosphere Pressure: 1.72 cm H₂O

Do not expose the mask to temperatures above 122°F (50°C).

BEFORE USE

- Inspect mask daily or prior to each use. Replace the mask if vent holes are blocked, if parts are damaged, or if gel is exposed due to tears or punctures.
- Check the anti-asphyxia valve. With the positive pressure device turned off, verify that the valve flapper is positioned so that room air can flow through the large opening in the valve (Figure 7a). Turn on the system; the flapper should close, and system air should flow through the mask (Figure 7b). If the flapper does not close or does not function properly, replace the mask. Do not block the opening on the anti-asphyxia valve or the vent holes. Ensure that the valve is not blocked with secretions and that the flapper is dry.

MOJO® 2 MASK FITTING INSTRUCTIONS

1. Verify the mask is the right size. When using the sizing guide ensure mouth is slightly open. Line the top of the guide to the bridge of the nose and select the smallest size that does not restrict the nose or mouth. The appropriate size may be smaller than expected due to the design.
2. Lay headgear flat with light gray side facing up with the shorter straps at the top. Place mask in the center of the headgear, and thread all 4 straps (Figure 1).
3. Unhook one or both magnetic clips from the mask. Seat the bottom cushion of the mask at the crease of the chin with your mouth slightly open and tip the mask up to contact the nose (Figure 2). Turn the forehead adjustment screw so the forehead pad lightly touches your forehead (Figure 6). The mask may sit lower on the nose than expected, depending on facial features.
4. Slip the headgear over the head, hook the magnetic clips in place. (Figure 3). The headgear should be flat across the back of the head.
5. Adjust lower headgear straps by using the hook and loop fastener tabs (Figure 4).
6. Adjust upper headgear straps: detach the hook and loop fastener tabs and pull straps forward to center mask over face. Ensure headgear is centered by aligning the two straps. Then pull the straps back, adjust strap tension and affix tabs. (Figure 5)
7. Turn the forehead pad adjustment screw to adjust tension until comfortable. If necessary, readjust top straps for a comfortable fit.
8. If the mask is uncomfortable repeat steps 5-8.

NOTE:

- I. If breathing through your nose or mouth is still constricted, contact your healthcare professional.
 - II. The mask should rest comfortably against the face. If there is excessive pressure on the bridge of the nose, loosen the top straps, slowly turn the adjustment screw counter-clockwise to relieve the pressure on the bridge (Figure 6). Then readjust the top straps and affix the tabs.
9. Attach the mask elbow to the device (CPAP, bi-level, ventilator) tubing; turn the device on and verify airflow through the mask.
 10. Gently press edges of the mask to change mask shape, create a comfortable fit, and reduce leaks.
 11. If leaks are noticed around the chin or the cheeks, slightly tighten the bottom straps.
 12. If leaks are noticed at the bridge of the nose, slowly turn the adjustment screw clockwise to slightly move the pad away from the forehead (Figure 6). Then tighten the top straps slightly until leaks are eliminated. This may be done in conjunction with molding the mask in the bridge area as described in Step 10 above.

NOTE:

- I. DO NOT overtighten any straps. Overtightening can worsen leaks. For optimal seal and comfort, mask should be snug, but not tight.

QUICK RELEASE OF MASK AND HEADGEAR

- Unhook magnetic clips and lift away from face. Remove mask and headgear assembly.

Mojo® is a trademark of Sleepnet Corporation.
Visit our website at www.sleepnetmasks.com.

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Medical Product Services GmbH