***Phantom<sup>®</sup>2 Nasal Vented Mask***

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The *Phantom*<sup>®</sup>2 Nasal Vented Mask is intended to be used with positive airway pressure devices, such as CPAP or bi-level, operating at or above 3 cm H<sub>2</sub>O. The mask is to be used on adult patients (>30kg) for whom positive airway pressure therapy has been prescribed. The mask is intended for single-patient multi-use in the home, hospital or institutional environment.

#### **NOTES**

- The mask system does not contain DEHP or phthalates.
- This product is not manufactured with natural rubber latex.
- If the patient has ANY reaction to any part of the mask system, discontinue use.
- An exhalation port is built into this mask so a separate exhalation port is not required.

#### **CONTRAINDICATIONS**

Do not use this mask if you or anyone (example: household members, bed partners, caregivers, etc.) in close physical contact with your mask has an active medical implant or metallic implant that will interact with magnets. Implant examples include, but are not limited to, pacemakers, implantable cardioverter defibrillators (ICD), neurostimulators, aneurysm clips, metallic stents, ocular implants, insulin/infusion pumps, cerebral spinal fluid (CSF) shunts, embolic coils, metallic splinter, implants to restore hearing or balance with implanted magnets (such as cochlear implants), flow disruption devices, contact lenses with metal, dental implants, metallic cranial plates, screws, burr hole covers, bone substitute device, magnetic metallic implants/electrodes/valves placed in upper limbs, torso, or higher, etc. If you have any questions regarding the implant, consult your physician or the manufacturer of your implant.

#### **CAUTIONS**

- U.S. Federal Law restricts this device for sale by or on the order of a physician.
- Do not expose the *Phantom*<sup>®</sup>2 Nasal Vented Mask to temperatures above 122°F (50°C).
- This mask is intended for single patient use only. It may be cleaned and used repeatedly on the same person but must not be used on multiple persons. Do not sterilize or disinfect.

#### **WARNINGS**

- This mask is not suitable for providing life support ventilation.
- Do not use the *Phantom*<sup>®</sup>2 Nasal Vented Mask with baseline pressure less than 3 cm H<sub>2</sub>O.
- If oxygen is used with this CPAP mask, the oxygen flow must be turned off when the CPAP machine is not operating. Explanation of the warning: When the CPAP device is not in operation, and the oxygen flow is left on, oxygen delivered into the ventilator tubing may accumulate within the CPAP machine enclosure. Oxygen accumulated in the CPAP machine enclosure will create a risk of fire.
- At a fixed rate of supplemental oxygen flow, the inhaled oxygen concentration will vary, depending on the pressure settings, patient breathing pattern, mask selection, and the leak rate.
- No smoking or open flames such as candles may be used when oxygen is in use.
- To minimize the risk of vomiting during sleep, avoid eating or drinking for three (3) hours before using the mask.
- Any unusual skin irritation, chest discomfort, shortness of breath, gastric distension, abdominal pain, belching or flatulence from ingested air, or severe headache experienced by the patient during or immediately after use should be immediately reported to a healthcare professional.
- Using a mask may cause tooth, gum or jaw soreness or aggravate an existing dental condition. If symptoms occur, consult a physician or dentist.
- Consult a physician if any of the following symptoms occur while using the mask: drying of the eyes, eye pain, eye infections, or blurred vision. Consult an ophthalmologist if symptoms persist.
- Do not overtighten the headgear straps. Signs of overtightening include: excessive redness, sores, or bulging skin around the edges of the mask. Loosen headgear straps to alleviate symptoms.
- At low CPAP pressures, the flow through the exhalation ports may be inadequate to clear all exhaled gas from the tubing. Some rebreathing may occur.
- This mask should be used with the positive airway pressure device recommended by a physician or respiratory therapist. The mask should not be used unless the positive airway pressure device is turned on and operating properly. Do not block or try to seal the exhalation port (vent holes). The vent holes allow a continuous flow of air out of the mask. When the positive airway pressure device is turned on and functioning properly, new air from the device flushes the exhaled air out through the mask's vent holes. When the positive airway pressure device is not operating, exhaled air may be rebreathed. Rebreathing of exhaled air can, in some circumstances, lead to suffocation. This warning applies to most models of CPAP machines and masks.

- This mask is not intended for patients without a spontaneous respiratory drive. This mask should not be used on patients who are uncooperative, obtunded, unresponsive, or unable to remove the mask. This mask may not be suitable for persons with the following conditions: impaired cardiac sphincter function, excessive reflux, impaired cough reflex, and hiatal hernia.
- Significant leak may occur without properly fitting the mask.
- Magnets are used in the mask and headgear clips with a field strength of 380mT. With the exception of the devices identified in the contraindication, ensure that the mask is kept at least 6 inches (approx.16 cm) away from any other medical implants or medical devices that can be impacted by the magnetic fields to avoid possible effects from localized magnetic fields. This applies to you or anyone in close physical contact with your mask.

### **GENERAL INFORMATION**

Dispose of the mask and accessories in accordance with local and regional guidelines.

This mask has a non-transferable three month limited warranty for manufacturer’s defects from the date of purchase by the initial consumer (or maximum service life, if shorter). If a mask used under normal conditions should fail, Sleepnet will replace the mask. For additional information about the Sleepnet Warranty, please visit <http://www.sleepnetmasks.com/>.

Per EU MDR 2017/745, please report serious incidents to manufacturer and EU Member State competent authority.

**NOTE:** The mask cushion has a service life of the 1 month. The swivel frame/tube set has a service life of 3 months. The TruFit Headgear and magnetic clips have a service life of 6 months.

### **CLEANING AND MAINTENANCE**

Cleaning is recommended on a daily basis. Refer to disassembly instructions. Wash the mask and components in warm water using a mild soap such as Ivory®. Rinse the mask thoroughly after washing. Visually verify that the vent holes are clear and allow the mask to air dry. Washing your face prior to fitting the mask will remove excess facial oils and help extend the life of the mask.

**NOTE:** The gel cushion, although silky to touch, will not tear, leak or detach from the mask under normal use and following recommended cleaning and maintenance protocols. However, the gel cushion is not indestructible. It will tear, cut and rip if subjected to mishandling or abuse. Do not clean in a dishwashing machine. The warranty does not cover these types of product failures. Therefore, please handle your new *Phantom*®2 Nasal Vented Mask with care.

### **TECHNICAL INFORMATION**

Resistance to air flow at 50 L/min: 0.45 cm H<sub>2</sub>O

Resistance to air flow at 100 L/min: 2.08 cm H<sub>2</sub>O

DEAD SPACE VOLUME (mL approx.): Mask – 30.7ml Tube – 66ml

Operating Temperature: Do not expose the mask to temperatures above 122°F (50°C).

### **BEFORE USE**

- Inspect mask daily or prior to each use. Replace the mask if vent holes are blocked, if parts are damaged, or if gel is exposed due to tears or punctures.

### **PHANTOM®2 FITTING INSTRUCTIONS**

1. Take mask out of the bag. Adjust the base strap according to your general size by aligning the tab at the indicated size markings (S, M, L, XL). (Figure 1).
2. Detach one magnetic clip from the mask. (Figure 2).
3. After placing the mask on your head, reattach the magnetic clip. (Figure 3).
4. First, adjust the top strap to raise or lower the mask. (Figure 4).
5. Second, adjust the left/right straps to make sure the bottom of the mask is sealed comfortably between the bridge of the nose and just above the upper lip. (Figure 5).
6. Connect the mask to the possible airway pressure device tubing and turn the device on to the prescribed pressure. Turn head from side to side while lying on the bed to assure secure fit. (Figure 6).
7. If unacceptable leaks occur, try moving the mask higher or lower around the nose to find the best fit. (Figure 7).
8. Gently press edges of the mask to change mask shape, create a comfortable fit, and reduce leaks. (Figure 8).
9. If leaks persist, try moving the mask higher or lower around the nose to find the best fit. (Figure 9).

### **DISASSEMBLY INSTRUCTIONS**

1. Firmly hold the mask cushion in one hand to one side. (Figure 1).
2. Grasp the swivel frame by one of the magnets. (Figure 2).
3. Firmly pull to separate the Swivel Frame from the cushion in a peeling motion. (Figure 3).

## **ASSEMBLY INSTRUCTIONS**

1. Lay cushion flat on a clean surface and oriented as shown (vent holes facing left). (Figure 1).
2. Align and fully insert the bottom Mask Post into the swivel frame slot as shown. (Figure 2).
3. Align the opposite post with the other Swivel Frame slot. Place one finger in the area shown and place another finger on the opposite side of the slot. Press down firmly where indicated until the 2nd post is inserted into the Swivel Frame. (Figure 3).
4. Turn the cushion over. Starting from the bottom, press the seal ring under the Swivel Frame tab. Work your way around the seal ring continuing to press the seal ring so that it seats around the Swivel Frame. The outer surface of the cushion should be flush with Swivel Frame. (Figure 4).

*Phantom*<sup>®2</sup> is a trademark of Sleepnet Corporation.

Visit our website at [www.sleepnetmasks.com](http://www.sleepnetmasks.com).

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Medical Product Services GmbH