



CPAP & NIV

Mojo[®] Full Face Mask

A mask that fits your life

The Mojo[®] Full Face Mask is an adjustable, soft-shell mask that delivers a personalized fit, functionality, and comfort to those who need it.



Available Sizes:

S M L XL

Custom Fit Technology

Sleepnet's Custom Fit Technology is our proprietary technology made from a flexible, pliable material that responds to gentle manipulation and accommodates changes in weight and facial asymmetry.

AIR^ogel[®]

Sleepnet's AIR^ogel[®] is an advanced cushion technology that conforms to the contours of each individual face creating a secure and comfortable seal.

Available Elbow Options

- Vented
- Non-Vented



AIR^ogel[®] Forehead Cushion

Our forehead cushion uses AIR^ogel[®] technology to provide enhanced comfort for the wearer.

Adjustable Tension Wheel

With just a turn of the wheel users can increase or decrease the tension of the mask, helping to redistribute pressure away from the nose bridge.

Easy-Release Magnetic Clips

Our magnetic clips are designed to simplify the process of attaching and detaching the mask from the headgear.

360[°] Rotating Elbow

Allows user to change sleep positions, move their head, or adjust the mask without discomfort.

What makes our Mojo® different?

Comfort - Sleepnet's AIR°gel® is a leading edge innovation developed and designed to enhance the comfort and effectiveness of our respiratory masks.

Personalized fit - Custom Fit Technology allows the user to easily shape the mask by gently bending and molding it along the sides, top, and bottom to provide a consistent and secure fit.

Adjustable - The Mojo® is designed with an adjustable tension wheel and AIR°gel® Forehead Cushion that are exclusive to the mask.

Easy on and off - our Easy-release magnetic clips are a user-friendly innovation designed to simplify the process of attaching and detaching the mask from the headgear.



Color coded elbow options for easy identification to meet your homecare and hospital needs.



Vented Mask with an anti-asphyxia valve for use with CPAP or bi-level devices.



Non-Vented Mask for use with dual limb circuits.

Mojo® Full Face Mask Specifications:

Product Number	Description	Size	Leak Port	Entrainment Valve	Mask Weight
50835	Mojo® Full Face Mask Vented with Headgear	S	✓	✓	144.0 g
50836	Mojo® Full Face Mask Vented with Headgear	M	✓	✓	148.0 g
50837	Mojo® Full Face Mask Vented with Headgear	L	✓	✓	156.0 g
50838	Mojo® Full Face Mask Vented with Headgear	XL	✓	✓	161.5 g
50920	Mojo® Full Face Mask Non-Vented with Headgear	S	N/A	N/A	144.0 g
50921	Mojo® Full Face Mask Non-Vented with Headgear	M	N/A	N/A	148.0 g
50922	Mojo® Full Face Mask Non-Vented with Headgear	L	N/A	N/A	156.0 g
50923	Mojo® Full Face Mask Non-Vented with Headgear	XL	N/A	N/A	161.5 g
55045	Mojo® Full Face Headgear only	One Size	N/A	N/A	N/A



Also by Sleepnet®
Aura Nasal

Contact us today:

- ☎ 1-800-742-3646
- ✉ info@sleepnetcorp.com
- 🌐 www.sleepnetmasks.com



The following product is protected by US patents 10,974,007 and 7,938,117.



Does not contain natural rubber latex

M90002 Version 0

