



CPAP

Aura[®] Nasal Mask

Entering a new era of comfort

The Aura[®] Nasal Mask is Sleepnet's first nasal mask featuring Custom Fit Technology and an Active Headgear Connector. Available in three sizes with a Touchless Spacebar, the Aura[®] will give you the comfort and rest you deserve.



Custom Fit Technology


Sleepnet's Custom Fit Technology is our proprietary, one-size-fits-all technology made from a flexible, pliable material that responds to gentle manipulation and accommodates changes in weight and facial asymmetry.

Quick-Release

Our quick release masks feature a slot on the side that allows the user to easily put on or take off the mask and headgear.

360° Rotating Elbow

Designed to provide flexibility and comfort, the 360° rotating elbow allows the user to change sleep positions, move their head, or adjust the mask without discomfort.



Touchless Spacebar

Adding comfort with its contactless design, our headgear attaches to the touchless spacebar and provides additional stability by keeping the mask seal in place when the user moves.

Available Sizes:



AIR[°]gel[®]

Sleepnet's AIR[°]gel[®] is an advanced cushion technology that conforms to the contours of each individual face creating a secure and comfortable seal.

Active Headgear Connector

Ideal for active sleepers, our active headgear connector is designed to move and adjust while maintaining a tight yet comfortable seal and preventing air leaks.



What makes our Aura® Mask different?

Comfort - One of the most significant advantages of AIR°gel® technology is the unparalleled comfort it provides. The gel molds to the user's face, reducing pressure points and maximizing comfort.

Better Seal - the Aura® mask is designed to accommodate most facial features and provide superior stability while maintaining a comfortable seal.

Ease of Movement - Our active headgear connector is designed to move and adjust as the wearer moves, allowing the mask to maintain a seal and prevent air leaks. This can be especially helpful for active sleepers, allowing for ease of movement during sleep.

Easy on and off - our quick release feature is a user-friendly innovation designed to simplify the process of putting on and taking off the mask and headgear.



AIR°gel® Technology:

Sleepnet's AIR°gel® is a leading edge innovation developed and designed to enhance the comfort and effectiveness of our respiratory masks. This advanced cushion technology features medical grade silicone gel that conforms to the contours of each individual face creating a secure and comfortable seal that ensures a snug fit while minimizing air leakage around the edges. Our AIR°gel® technology is soft and flexible and helps to reduce irritation, making it suitable for extended wear.




Aura® Nasal Mask with EZFit Headgear Specifications:

Product Number	Description	Leak Port	Cushion Size	Mask Weight with Cushion
50176	Aura® Vented Nasal Mask	✓	S	61.5 g
50177	Aura® Vented Nasal Mask	✓	M	62.5 g
50158	Aura® Vented Nasal Mask	✓	L	71.5 g
55051	EZFit Headgear only	N/A	N/A	N/A



Also by Sleepnet®
Mojo® Full Face

Contact us today:

-  1-800-742-3646
-  info@sleepnetcorp.com
-  www.sleepnetmasks.com



The following product is protected by US patents 10,974,007, 9,199,054, and other patents pending



FDA
Cleared

Does not contain
natural rubber latex

M90017 Version 0



Sleepnet