





## What makes our Ascend® Mask different?

**Comfort** - One of the most significant advantages of AIR°gel® technology is the unparalleled comfort it provides. The gel molds to the user's face, reducing pressure points and maximizing comfort.

**Better Seal** - the Ascend mask is designed to accommodate most facial features and provide superior stability while maintaining a comfortable seal.

**Ease of Movement** - Our active headgear connector is designed to move and adjust as the wearer moves, allowing the mask to maintain a seal and prevent air leaks. This can be especially helpful for active sleepers, allowing for ease of movement during sleep.

**Available for resupply** - Our replaceable cushions can be removed and cleaned, helping to maintain better hygiene and mask performance.

### Easy-to-Replace Cushion S M L




**Assembly:** Place the top alignment pin in the notch. Push the bottom until all pins are secure.



**Disassembly:** Squeeze the sides of the cushion frame to release from mask.

### Ascend® Full Face Mask with AFit Headgear Specifications:

Product Number	Description	Cushion Size	Mask Weight with Cushion
50825	Ascend® Full Face Mask system and headgear	S,M,L	N/A
50834	Ascend® Full Face Mask system with attached headgear	S,M,L	N/A
50801	Ascend® Full Face Mask with attached headgear	S	94.5 g
50802	Ascend® Full Face Mask with attached headgear	M	103.4 g
50803	Ascend® Full Face Mask with attached headgear	L	107.8 g
50813	Ascend® Full Face Mask without headgear attached	S	94.5 g
50814	Ascend® Full Face Mask without headgear attached	M	103.4 g
50815	Ascend® Full Face Mask without headgear attached	L	107.8 g
50207	AIR°gel cushion	S	NA
50208	AIR°gel cushion	M	NA
50209	AIR°gel cushion	L	NA
55058	AFit headgear only	N/A	N/A



Also by Sleepnet®  
**Ascend Nasal**

#### Contact us today:

☎ 1-800-742-3646  
✉ info@sleepnetcorp.com  
🌐 www.sleepnetmasks.com



The following product is protected by US patents 10,974,007, 9,199,054, and other patents pending



FDA  
Cleared

Does not contain natural rubber latex

M90056 Version 0



Sleepnet