



NIV

Veraseal[®]3 Full Face Mask

Acute care comfort

The Veraseal[®]3 disposable full face mask was designed for an acute care environment with a focus on patient comfort. Easily adjustable to the patient's face, the Veraseal[®]3 provides stability and ease of movement.

AIR[°]gel[®]

Sleepnet's AIR[°]gel[®] is an advanced cushion technology that conforms to the contours of each individual face creating a secure and comfortable seal.

Adjustment Dial

By turning the dial, users can increase or decrease the pressure exerted by the pad against their forehead to help stabilize the mask and redistribute pressure away from the bridge of the nose.

AIR[°]gel[®] Forehead Cushion

Our forehead cushion uses AIR[°]gel[®] technology to provide enhanced comfort for the wearer.

Quick-Release

Our quick release masks feature a slot on the side that allows the user to easily put on or take off the mask and headgear.

360° Rotating Elbow

Allows user to change sleep positions, move their head, or adjust the mask without discomfort.

Oxygen Port

Allows easy delivery of supplemental oxygen.

Available Sizes:

S M L XL



What makes our Veraseal®3 different?

Comfort - One of the most significant advantages of AIR°gel® technology is the unparalleled comfort it provides. The gel molds to the user's face, reducing pressure points and maximizing comfort.

Easy on and off - the Veraseal 3 provides a quick release slot that allows providers to easily put on or take off the mask and headgear.

Adjustable - The Veraseal®3 is designed with an adjustable dial that is exclusive to the mask.

Reduces skin breakdown - The Veraseal®3 is easily adjustable and built with Sleepnet's AIR°gel® cushion, which provides a comfortable seal while minimizing irritation.

Available in four sizes - Small, Medium, Large, and Extra Large.



Color coded elbow options for easy identification to meet your hospital needs.



Vented with an anti-asphyxia valve for use with CPAP or bi-level devices.



AAV Non-Vented with an anti-asphyxia valve for use with single limb circuits.



Non-Vented for use with dual limb circuits.

Veraseal®3 Full Face Mask Specifications:

Product Number	Description	Size	Leak Port	Entrainment Valve	Mask Weight
50843	Veraseal®3 Vented Full Face Mask	S	✓	✓	109.2 g
50844	Veraseal®3 Vented Full Face Mask	M	✓	✓	114.3 g
50845	Veraseal®3 Vented Full Face Mask	L	✓	✓	119.4 g
50846	Veraseal®3 Vented Full Face Mask	XL	✓	✓	120.1 g
50949	Veraseal®3 Non-Vented Full Face Mask	S	N/A	N/A	108.6 g
50950	Veraseal®3 Non-Vented Full Face Mask	M	N/A	N/A	113.7 g
50951	Veraseal®3 Non-Vented Full Face Mask	L	N/A	N/A	118.8 g
50952	Veraseal®3 Non-Vented Full Face Mask	XL	N/A	N/A	119.5 g
50953	Veraseal®3 AAV Non-Vented Full Face Mask	S	N/A	✓	109.1 g
50954	Veraseal®3 AAV Non-Vented Full Face Mask	M	N/A	✓	114.2 g
50955	Veraseal®3 AAV Non-Vented Full Face Mask	L	N/A	✓	119.3 g
50956	Veraseal®3 AAV Non-Vented Full Face Mask	XL	N/A	✓	120.0 g



Also by Sleepnet®
Innova® Nasal

Contact us today:

- 1-800-742-3646
- info@sleepnetcorp.com
- www.sleepnetmasks.com



The following product is protected by US patents 10,974,007, 9,199,054, and D892,304.



FDA
Cleared

Does not contain natural rubber latex

M90082-01 Version 0

