



CPAP & NIV

# iQ<sup>®</sup> Blue Nasal Mask

## Smart technology for better therapy

A smart choice for a custom mask, our iQ<sup>®</sup> Blue Nasal mask adjusts to the contours of your face, providing optimal fit with maximum comfort.


### Custom Fit Technology

Sleepnet's Custom Fit Technology is our proprietary technology made from a flexible, pliable material that responds to gentle manipulation and accommodates changes in weight and facial asymmetry.



### AIR<sup>°</sup>gel<sup>®</sup>

Sleepnet's AIR<sup>°</sup>gel<sup>®</sup> is an advanced cushion technology that conforms to the contours of each individual face creating a secure and comfortable seal.



### Quick-Release

Our quick release masks feature a slot on the side that allows the user to easily put on or take off the mask and headgear.

### Available Elbow Options

- Vented
- Non-Vented

### 360° Rotating Elbow

Designed to provide flexibility and comfort, the 360° rotating elbow allows the user to change sleep positions, move their head, or adjust the mask without discomfort.

## What makes our iQ® Blue mask different

**Comfort** - Sleepnet's AIR°gel® is a leading edge innovation developed and designed to enhance the comfort and effectiveness of our respiratory masks.


**Flexible** - the iQ Blue Nasal mask can be shaped to each individual face, providing flexibility and extended visibility.

**Personalized fit** - Custom Fit Technology allows the user to easily shape the mask by gently bending and molding it along the sides, top, and bottom to provide a consistent and secure fit.

**Clear Field of Vision** - The iQ Blue mask is designed to offer a Clear Field of Vision giving the user the ability to watch television, read, or wear glasses while the mask stays comfortably in place.



## The iQ® Blue offers two headgear choices:



### 3-Point Headgear

With a minimalist design, the 3-Point headgear offers three adjustment points for a comfortable fit.



### StableFit Headgear

The StableFit Headgear offers five adjustment points that allow for a more secure fit.

## iQ® Blue Mask Specifications:

Product Number	Description	Leak Port	Headgear Style	Mask Size	Mask Weight
50647	iQ® Blue (Vented) Moldable Nasal Mask w/ Headgear	✓	3-Point	One Size	63.8 g
50655	iQ® Blue (Vented) Moldable Nasal Mask w/ Headgear	✓	Stablefit	One Size	63.8 g
50660	iQ® Blue (Non-Vented) Moldable Nasal Mask w/Headgear	N/A	3-Point	One Size	63.9 g



Also by Sleepnet®  
**iQ<sup>2</sup>**

### Contact us today:

☎ 1-800-742-3646  
✉ info@sleepnetcorp.com  
🌐 www.sleepnetmasks.com



The following product is protected by US patents 10,974,007, 7,938,117 and other patents pending



Does not contain natural rubber latex

M90004 Version 1

