



CPAP & NIV

iQ[®] Blue Nasal Mask

Smart technology for better therapy

A smart choice for a custom mask, our iQ[®] Blue Nasal mask adjusts to the contours of your face, providing optimal fit with maximum comfort.


Custom Fit Technology

Sleepnet's Custom Fit Technology is our proprietary technology made from a flexible, pliable material that responds to gentle manipulation and accommodates changes in weight and facial asymmetry.



AIR[°]gel[®]

Sleepnet's AIR[°]gel[®] is an advanced cushion technology that conforms to the contours of each individual face creating a secure and comfortable seal.



Quick-Release

Our quick release masks feature a slot on the side that allows the user to easily put on or take off the mask and headgear.

360° Rotating Elbow

Designed to provide flexibility and comfort, the 360° rotating elbow allows the user to change sleep positions, move their head, or adjust the mask without discomfort.

What makes our iQ® Blue mask different

Custom Fit Technology - Our proprietary, one-size-fits-all technology that responds to gentle manipulation and accommodates changes in weight and facial asymmetry.

AIR°gel® Cushion Technology - An advanced cushion technology that conforms to the contours of each individual face creating a secure and comfortable seal.

Clearest Field of Vision Possible - No forehead pad required

Advanced Cushion Technology - Our patented AIR°gel technology with an ergonomic design and luxurious feel reduces pressure while maintaining an effective seal.



The iQ® Blue offers two headgear choices:



3-Point Headgear

With a minimalist design, the 3-Point headgear offers three adjustment points for a comfortable fit.




StableFit Headgear

The StableFit Headgear offers five adjustment points that allow for a more secure fit.

iQ® Blue Mask Specifications:

Product Number	Description	Leak Port	Headgear Style	Mask Size	Mask Weight
50647	iQ® Blue (Vented) Moldable Nasal Mask w/ Headgear	✓	3-Point	One Size	63.8 g
50655	iQ® Blue (Vented) Moldable Nasal Mask w/ Headgear	✓	Stablefit	One Size	63.8 g



Also by Sleepnet®
iQ2

Contact us today:

☎ 1-800-742-3646
✉ info@sleepnetcorp.com
🌐 www.sleepnetmasks.com



The following product is protected by US patents 11,998,694, 10,974,007, 7,938,117 and other patents pending



FDA
Cleared

Does not contain natural rubber latex

M90001 Version 2

